



EASYCAP

EEG Recording Caps and Related Products

EASYCAP GmbH
Am Anger 5
82237 Würthsee
GERMANY

Tel +49 8153 88702-00
www.easycap.de
office@easycap.de


32Ch Standard Cap for actiCAP

 Green Holders: Label 1-32, Hardware Channel 1 – 32

 Blue holder: Label & hard-wired Ref

 Black holder: Label & hard-wired Gnd

Electrode Names and Number Labels



Details

Ordering Information

For ordering please give **Article Number**, **Cap Cut**, and **Size**

(Example: CAC-32, C-Cut, 58):

- Article Number: **CAC-32**
- Cap Cut: **C-Cut** or **A-Cut**
- Size (given in cm head circumference):
 - Adult caps: **54, 56, 58, 60, 62, 64** (average male: 58, average female: 56)
 - Children caps: **50** (3-4 years), **52** (5-10 years), **54** (11-14 years)
 - Infant caps: **34, 36** (newborn), **38, 40** (3 months), **42, 44** (7 month), **46, 48** (2 years)
 - Preterm caps: **32, 34**

The catalogue-number comprises the cap as described. For further information about accessories or consumables, please visit our website or contact our local distributor.

Cap

White Subtemporal Cap with integrated chin belt,

Sizes 52 – 64 made from High Precision Fabric, Sizes 50 and smaller made from High Comfort Fabric

Options: C-Cut or A-Cut, Size. For further variation please contact us.

The labelled holders are mounted according to the figure (1-32 transparent, Ref blue, Gnd black). The positions include the 19 sites of the "International 10/20-System". The core positions are extended by further positions of the 10%-System to increase the topographical resolution. Positions TP9/TP10 approximately meet left and right mastoids.

Each cap comes with 3 additional holders (apricot empty label) for use with double-sided adhesive rings to place electrodes on bare skin. The grip of the washer is better, when both, skin and holder, are degreased with alcohol prior to application.

Accessories

The cap can be fixed either by the integrated chin belt alone, or additionally / alternatively by a chest belt (not included). Normally the decision between these two methods can be made according to what is more comfortable for the subject. Under difficult conditions the chest belt will hold the cap in place a little better.

Service

In case of any further questions please contact our local distributor or support@easycap.de . Please have the serial number of your cap at hand (syntax Axxxxx, located on the backside label of the cap).