


Wireless Trigger

Wireless Trigger transmission made easy



- ✓ Wirelessly transmits event markers to your EEG amplifier
- ✓ Guarantees fast and consistent trigger transmission
- ✓ Allows sending triggers to multiple receivers simultaneously
- ✓ Encodes your stimulus events with 8-bit triggers
- ✓ No software configuration necessary



Wireless Trigger: Wireless recording of event markers in EEG data

- ✓ Fully mobile solution for recording event markers in EEG data
- ✓ Transmit 8-bit trigger information to any compatible amplifier
- ✓ The transmitter is plugged into the parallel port of your PC or into the TriggerBox output. Use the STE compact as a receiver or connect Wireless Trigger Receivers to 9-pin DSub input ports
- ✓ This solution provides the fastest and most consistent wireless trigger timing
- ✓ With just 2.5 ms transmission delay and only one sample jitter, it has the most accurate timing that is needed for gold standard EEG data analysis
- ✓ It is the ideal companion for your mobile EEG recordings with our portable LiveAmp amplifiers

